



E-Walk

ENG

2018/2019



Easy and
convenient

MADE BY DECON
SINCE
Creating Freedom
1989
POWERED BY YOU

balance

24 < 46 >

← LADDUTTAG

→ decon wheel


E-Walk

E-WALK

Simple drive unit

WHAT DO I GET WHEN I ORDER E-WALK?

- Drive unit
- Battery
- Battery charger
- Operating unit
- User manual
- Two-year warranty



E-Walk is an attendant-controlled drive unit that provides additional traction when driving the wheelchair uphill and braking when coming downhill.

When E-Walk is not needed, the motor can be removed in one simple action. Control unit, with comfortable handle (not heavy, easy to hold) mounted directly on the wheelchair handle. E-walk can easily be used for different chairs at the same time. Up to 25 km distance with one battery charge.

E-walk is just as suited to indoor and outdoor use.

TECHNICAL DATA

Motor	24V 200 W
Battery type	LiFePo4 24V15Ah/lead 2x12v 12Ah
Total weight	15.4 kg (without lead battery)
Maximum user weight	Maximum 150 kg
Suitable for chair width	38–60 cm
Drive wheel	200x50, Solid
Front wheels	75x22, Solid
Forward speed	0–6 km/h
Reverse speed	0–2.6 km/h
Maximum range	Up to 22 km in accordance with ISO 7176-4
Fuse	25 A (flat)

Battery charger

Capacity	100 to 240V AC, 50/60Hz
Charging time	3 hours LiFePO4/5–6 hours lead
CE mark	In accordance with EU Directive 93/42/EEC
User instructions	See user manual

ITEM NUMBER E-WALK

LEAD BATTERY	MCG4050
LIFEPO4 BATTERIES	MCG4750

LiFePO4 batteries (Lithium iron phosphate) are rechargeable batteries which have been technologically developed from lithium ion batteries. They are lightweight, have a high power output and slow self-discharge rate. LiFePO4 batteries can handle up to 1,600 cycles.



Please feel free to contact us
on +46 (0)345 408 80
or at info@decon.se

ACCESSORIES

Moving control panels

380060

If there is no space due to handles for other functions, simply move the throttle grip.



It can be moved to either the left or right side.

Battery

Lead battery

MCG2001

Battery charger for lead battery (*Two are required per machine*)

MCG2002

LiFePO4 battery

MCG2015-1

Battery charger for LiFePO4 batteries

MCG2016

Throttle control

Fixed cord

MCG2003

Quick connection cord

MCG2003-1



Converter

MED12V

Charge your battery in the car, boat, caravan or motor home
Output 230VAC 50Hz 150W



Seat

E-Walk is equipped for seat width 38–45 cm as standard.

If the wheelchair has a seat width wider than 45 cm, E-Walk can be supplemented with wider side pieces.

For seat widths 45–53 cm

350165B

The complete side piece, right and left, E-Walk, width 45–53

For seat width 53–60 cm

350165C

The complete side piece, right and left, E-Walk, width 53–60 cm



ADAPTERS

	ITEM NO.	PRICE CAT.
Alu Rehab		
Netti I/Netti II/Netti III	MCG1107	40
Bescomed		
Primus ML	MCG1149	10
Berollka Active		
Active Basic	MCG1140	30
Active Finess	MCG1113	40
B2	MCG1173	20
F1	MCG1150	30
Rookie	MCG1172	20
Bischoff & Bischoff		
Pyro Light XL	MCG1144	30
Pyro Start (22" and 24")	MCG1162	20
S-Eco 2	MCG1161	20
Dietz		
Primo Amico	MCG1139	20
Primo Amico (new)	MCG1118	20
Caneo E/s	MCG1163	20
Drive Medical		
Enigma	MCG1164	30
Etac		
Cross	MCG1106	40
Cross 5	MCG1171	40
Cross XL	MCG1125	40
M100TR	MCG1188	20
M100	MCG1187	
Next/Next Comfort	MCG1121	40
Transit	MCG1105	30
Twin	MCG1117	50
Handicare		
Emineo (Basic & One Tool Mounting)	MCG1114*	30
Exigo 20/Legend 2 Exigo	MCG1156	20
Ibis XA	MCG1119	30
Ibis XA 2013	MCG1178	20
Legend 2 (One Tool Mounting)	MCG1115	20
Legend 2 (Basic Mounting)	MCG1116	20
HD Rehab		
HD 500	MCG1109*	30
HD 500 (20" with 125 mm castors)	MCG1154	30
Humlan 500	MCG1109	30
HD 600	MCG1110	30
HD 650	MCG1111	30
Humlan 650	MCG1111	40
HD Balance (16" and 24")	MCG1189	40

*Low MODEL E-Walk (LiFePO4 batteries)



	ITEM NO.	PRICE CAT.
Invacare		
Action 3	MCG1165	20
Action 1Ng	MCG1193	20
Action 4 Ng	MCG1155	20
Adapt	MCG1122	20
Rea Assist	MCG1147	20
Rea Azalea	MCG1145	20
Rea Azalea Assist/Tall/Minor	MCG1145	20
Rea Azalea (new chassi)	MCG1192	
Rea Focus/Focus 150	MCG1123	30
Rea 703 Lx	MCG1133	30
Rea Spinx	MCG1130	30
Rea Spirea 2	MCG1128	10
Rea Spirea 4 Ng 16"	MCG1179	10
Küschall		
Compact	MCG1141	30
Life & Mobility		
Roux	MCG1151*	40
Meyra		
Active Line	MCG1157	10
Eurochair Avanti Model 1.736	MCG1166	10
Eurochair Basic Model 1.751	MCG1170	20
Eurochair Vario Model 1.750	MCG1170	20
Eurochair Model 1.850	MCG1143	30
Format Model 3.940	MCG1185	30
Format 2 Model 2.950	MCG1167	30
X3 Model 4.352	MCG1158	30
Solero Light Model 9.072	MCG1159	30
Minos		
Global	MCG1168	40
Irati alu	MCG1152	10
Plena	MCG1169	10
Otto Bock		
Innov'Effect	MCG1148	20
Motus cs/Motus cv	MCG1182	40
Start m2 effect	MCG1148	20
Start m5 comfort	MCG1184	40
Sunrise Medical		
Breezy Basix	MCG1137	20
Breezy Entrée	MCG1138	20
Breezy Helix	MCG1124	30
Breezy Parix	MCG1186	20
Breezy Rubix	MCG1124	30
Quickie Neon	MCG1120	30
Sopur Easy 160l	MCG1142	30
Sopur Easy 200/300	MCG1112	30
Sopur Starlight	MCG1135	20
Sopur Youngster	MCG1136	30
Van Os Medical		
Excel G-Eco	MCG11694	30
Vermeiren		
Jazz	MCG1160	20

**DO YOU HAVE ANY QUESTIONS
OR WANT TO ORDER?**

+46 (0) 345 408 80 – We'll be happy to help!

Checklist when ordering

- For standard equipment see "What do I get when I order an E-Walk"
- Ordering the necessary adapter
- You need to order one adapter kit for one E-Walk
- If you cannot find the necessary adapter in this adapter list, see our website for updated information. If this does not appear here either, please contact us
- Check your order confirmation, if anything is incorrect or unclear, please contact us immediately
- When you receive your products, please check the packaging for damage. In such case, contact us before you sign for receipt of the products



 **decon**

Södra Ekeryd 119 314 93 Hyltebruk

Telephone: +46 (0)345 40880 Fax: +46 (0)345 40895 E-mail: info@decon.se www.decon.se

Decon® is not responsible for any changes made after printing.

© All rights are reserved for Decon®. With reservation for printing and typographical errors, discontinued products and misleading images.