



# O'Neill Healthcare

Enhancing Quality of Life, Since 1999



## State of the Art Lifting Chair

Liftup Raizer is a battery operated emergency lifting device that helps a fallen person get from the floor to a sitting position in a safe and timely manner. Liftup Raizer is easy to operate by only one assistant and does not require any physical effort, besides a supportive hand

### Equal Opportunities

Every day, thousands of elderly and physically impaired people are falling without being able to get up without assistance. The Raizer from Liftup can be easily assembled and operated by only one assistant. The fallen person can be raised in a safe and effortless manner within minutes.

### Easy Operation

The assistant can safely operate the Raizer with minimal physical effort while at the same time observing and caring for the fallen person by, for example, supporting their head throughout the transfer process. The Raizer is an optimal solution for personnel in home care as well as for ambulance services and all personnel working with handling and transfer of individuals in general. The Raizer is light and easy to transport. It is carried disassembled to the place where it is to be used. The Raizer is easily assembled and placed around the

fallen individual with minimal effort. The Raizer is battery operated and has operation buttons on the side of the chair – and it can be supplied with a remote control as supplementary equipment.

### Safety above all

For a fallen individual it is of utmost importance that there is a feeling of security while the transfer to standing or sitting position is taking place. At all times during the transfer with Raizer the individual has contact with the floor as well as physical contact and eye contact with the assistant. For the assistant a transfer can be done without any strain to his or her back. The Raizer can transfer individuals weighing up to 150 kgs. The transfer movement is smooth and protective and as an extra comfort the Raizer is equipped with a seat belt.

### BENEFITS

- Mobile chair for transfer from floor to a sitting or standing position.
- Minimal physical effort for the assistant
- Both the fallen person and the assistant can feel comfortable and safe in the process
- Easy to transport in 2 light bags
- Assembly prior to operation: max 3 minutes




# O'Neill Healthcare

Enhancing Quality of Life, Since 1999



## Technical Specifications

- Lifting capacity: Max. 150 kg
- Lifting time: 20-30 secs.
- Weight:  
Seat 9 kg. + legs and back rest 4 kg  
– total 13 kg
- Charger (2 pcs.):  
Soniel 12V. 100-240V and  
12 V car adapter
- Number of lifts on a full charge:  
40 with a max. load Approx. 100 with an  
average load
- Charging time: Max. 6 hours
- Charging of empty battery:  
10–15 min = 1 transfer
- Easy to clean
- CE-mark
- Patent pending



## Extras

- Remote control
- Hygiene cover
- Trolley
- Stair wheels
- Bike bracket